



CIRCULATION ULTRA

circulation.co.za

Included in this standard package:

- Device x 1pc
- Medium pad x 4pcs
- Small pad x 2pcs
- Large pad*2pcs
- Electric wire x 3pcs
- (when use 4 line cable, make sure to use at least 1pc white(positive) +1pc blue line (negative) at the time.
- USB cable x 1set
- Acupuncture point chart x 1pc
- Instruction manual x 1pc
- Collection board for pad and wire x 1pc
- Touch pen x 1pc
- Bag x 1pc


Features

- 4" touch screen
- 2 output with 12 modes
- Dual Channel
- 20 level intensity
- 20 minutes default timer, can adjustable from 10-60 minutes
- Rechargeable lithium battery
- Battery indicator to know when need to recharge battery
- USB Charger 5v




1

Illustrations of the device & accessories: (Model: AS2015)



- 1.ON/OFF button
- 2.Battery indicator
- 3.Intensity(+)/Button
- 4.Intensity(-)/Button
- 5.LCD screen
- 6.Time setting button
- 7.Pause button
- 8.Two output connector
- 9.USB port
10. A/B shift button
- 11.Safety lock icon

2



- 12.Pads-sticking board
- 13.USB cable
- 14.Pads Wires
- 15.Massage Pads

Modes:

- There are twelve modes to choose from,there are:
- 1.Tuina
 - 2.Acupuncture
 - 3.Tapping
 - 4.Guasha
 - 5.Cupping
 - 6.Shiatsu
 - 7.Pressure
 - 8.Back Pain
 - 9.Deep Tissue Massage
 - 10.Body Building
 - 11.Weight Loss
 - 12.Combination

3

Mode	Name	Description	Function
1	Tuina	Simulating Tuina is a hands-on body treatment which opens the flow of Qi through the meridians.	It's especially effective for joint pain,sciatica,muscle spasms,and pain in the back,neck,and shoulders.It also helps chronic conditions such as insomnia,constipation,headaches,and the tension associated with stress.
2	Acupuncture	Simulating acupuncture provides the same sensation without needles,by utilizing electronic pulses to achieve the sensation of differing needle depths.	By stimulating specific acupuncture points,it corrects imbalances in the flow of qi through channels known as meridians.Please read the acupuncture chart.
3	Tapping	Simulating tapping is intended to produce a hard thumping sensation that massage your muscles,which can give you an invigorated and revitalized feeling.	It's especially effective for muscle relaxing.
4	Guasha	Simulating Guasha with a low frequency setting provides the user with a scratching and kneading sensation.Gua sha releases unhealthy elements from injured areas and stimulates blood flow and healing.	Guasha is effective in acute and chronic internal organ disorders.Patients may see immediate relief that patient feel from pain,stiffness,fever,chill,cough,wheeze,nausea and vomiting etc.
5	Cupping	Der simulierende Shiatsu ist eine alternative Medizin, bei welcher Druck von Fingern und Händen ausgeübt wird und der Körper so massiert wird.	Cupping is used to relieve back and neck pains,stiff muscles,anxiety,fatigue,migraines,rheumatism,and even cellulite etc.
6	Shiatsu	Simulating Shiatsu,it is a type of alternative medicine consisting of finger and palm pressure,stretches,and other massage techniques.	Shiatsu can help facilitate healing from many minor ailments.From headaches and back pain to digestive,sinus and skin problems.
7	Pressure	It is the manipulation of superficial and deeper layers of muscle and connective tissue using various techniques.	Can enhance function,aid in the healing process,decrease muscle reflex activity,inhibit motor-neuron excitability,promote relaxation and well-being,and as a recreational activity.
8	Back Pain	Use low frequency vibrations to promote circulation of blood and relieve pain.	Specifically use for lower back.
9	Deep Tissue Massage	Deep tissue massage therapy is similar to Swedish massage,but the deeper pressure is beneficial in releasing chronic muscle tension.	Deep tissue massage can provide relief from pain and stiffness associated with sports injuries,whiplash,osteoarthritis,muscle spasms,fibromyalgia and strain injuries such as carpal tunnel syndrome.
10	Body Building	Applying electrical current to muscles causes muscles to contract.Stimulating muscles repeatedly with electricity may eventually result in muscles that are strengthened and toned to some extent.	Electrical muscle stimulation helps muscle reduction,relaxation of muscle spasms,increase range of motion,injury recovery and post surgery recovery.
11	Weight Loss	Combination of acupuncture and other waveforms.	It helps lose weight or tone,strengthen and firm muscles at abs,back,thighs and arms etc.
12	Combination	Auto combination of several modes for relaxation.	Full body relaxation after workout:exercises&fitness or after a day of hard work.

4

Recommended massage points:




LOWER BACK
Modes 1,4,5,6,7,10,11



DELTOID
Modes 2,3,4
Attach one pad in front and one in back of the muscle.



ELBOW
Modes 1,2,3,4,6,7,9,11
Attach both pads on either side of the joint with the pain.



SCIATIC
Modes 1,2,3,4,5



HIP & THIGH
Modes 1,2,3,4,5,6,7,10
Attach both pads on either side of the area with pain.




HIP & THIGH
Modes 1,2,3,4,5,6,7
Attach both pads on either side of the area with pain.

6




KNEE
Modes 1,2,3,4,5,6,7,10
Attach both pads above the knee or above and below joint with pain.



ANKLE/FOOT
Modes 1,2,3,6,7,12
Attach pads per the illustration, on the left/ on the outside of your ankle/foot for pain.Attach the pads per illustration on the right / on the other side of your ankle/foot

7



CALF
Modes 1,2,4,5,6
Attach both pads on the calf/leg where you feel pain. (Do not place electrode pads simultaneously to the calves of both legs.)



ABS/SHOULDER
Modes 1,4,5,6,7,10,11

7

Included in this standard package:

1. Connect the pads to the wires by snapping them on, then remove The protective film. Note: adding a few drops of water on the pads will increase their longevity.
2. Place pads where the massage sensations are desired and plug into one of the ports on the device. For additional relief use the second port and repeat steps 1&2 for two additional pads. Make sure both pads are on the skin and are not overlapping. Make sure your skin is free from any dirt, oil, or lotion.
3. Turn the device on by pressing the on/off switch on the top of the unit. The unit will be now activated and start in Mode 1.
4. Touch the screen to choose modes. Then touch center button off the screen to control A/B channel. Remember the strength will default to the lowest setting with each mode selection.
5. Once you select a mode, gradually increase the intensity by pressing the "+" button; decrease it by pressing the "-" button on the side of main device.
6. Recommended application is 10-60 minutes per day. Increase 10 minutes each setting. If desired, the time remaining automatically counts down to zero and the unit will automatically turn off.

Features

- Do not apply the massage pads near the heart, on the head, above the neck, in the pubic region, over Scarred areas, on the throat or over the mouth.
- Avoid touching the pads when the unit is on.
- Pregnant women must avoid using this massager.
- Do not use this massager if you a pacemaker or other life

8

Included in this standard package:

1. Connect the pads to the wires by snapping them on, then remove The protective film. Note: adding a few drops of water on the pads will increase their longevity.
2. Place pads where the massage sensations are desired and plug into one of the ports on the device. For additional relief use the second port and repeat steps 1&2 for two additional pads. Make sure both pads are on the skin and are not overlapping. Make sure your skin is free from any dirt, oil, or lotion.
3. Turn the device on by pressing the on/off switch on the top of the unit. The unit will be now activated and start in Mode 1.
4. Touch the screen to choose modes. Then touch center button off the screen to control A/B channel. Remember the strength will default to the lowest setting with each mode selection.
5. Once you select a mode, gradually increase the intensity by pressing the "+" button; decrease it by pressing the "-" button on the side of main device.
6. Recommended application is 10-60 minutes per day. Increase 10 minutes each setting. If desired, the time remaining automatically counts down to zero and the unit will automatically turn off.

Easy steps to get started with your therapy

- Regular Application principles:
1. Find the exact pain point; this is the position for the pads to be placed. The acupuncture chart can also serve as a guide for pad placement specific to each ailment.
 2. Adjust the intensity gradually. It is more effective to reach the highest intensity level you can stand.
 3. Recommended application duration is 2-3 times per day for 40-60 minutes per application, over 10 days for one course.
 4. Pads can be placed in two ways:
 - Twin Mode-with the 2 pads next to each other on the same side of the body and in close proximity.
 - Opposed Mode-with the 2 pads on a similar spot, but on either sides of the body, i.e. the front and back shoulder.

Maintenance and care

1. When you are finished using the device, turn it off and take off the pads. Store the unit in a cool, dry place until next use.
2. After usage, please cover the pads with the protective film and put them back in their zip lock bags. Applying a few drops water before and after usage will increase the longevity of the pads.

9

3. Charge the unit at least every 3 month in event it is not used regularly. You can also use the unit with the wall charger or directly from your computer using the USB cable.

Troubleshooting

Problem	Cause	Solution
One pad feels stronger than the other.	This is normal. Different areas of your body will react differently.	Nothing needs to be done. Make sure the pads are moist and are making good contact.
During the massage, the skin feels a painful burning sensation or the stimulation becomes weakened.	Adhesive gel pads are not adhering firmly to the skin. Or the gel pads are too dry.	Apply a few drops of water to adhesive surface of each pad and make sure the pads are pressed firmly to the skin during application.
Unit is on, but no sensation is felt in the pads.	Both pads are not firmly on the skin. Or there is a loose connection. Or the intensity level needs to be increased.	Make sure both pads are firmly pressed to the skin. Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.
Adhesive gel pads do not stick to skin even after cleaning and moistening the gel pad.	Adhesive gel pads need to be replaced.	Replace the set of pads. Pads can be purchased from local distributor.
Unit does not turn on or the screen is dim.	Battery is low.	Recharge the battery.
Does not provide adequate pain relief.	Not using it long enough. Not lying down then massaging neck.	Use your massager 20-30 minutes at a time, 3-6 times per day. Lie on your back for neck support when you need relief for headaches or neck pain.

10

circulation.co.za